

prisma SOFT plus / prisma SMART plus prisma SOFT max / prisma SMART max

CPAP and APAP therapy: Flexible. Secure. Reliable.

With prisma SMART and prisma SOFT in the **plus** and max variants, we offer modern, reliable CPAP and APAP therapy featuring integrated connectivity, precise event detection and a high level of patient comfort.

The **plus** stands for integrated Bluetooth for data transmission to the prisma APP and thus for modern patient self-management. Equipped with an integrated modem, the **max** variants also offer the complete telemedical prisma CLOUD connection.

Our Software Solutions



prisma CLOUD

The telemedical connection for experts and patients. Information updated daily on therapy, fully automated, from anywhere.

• prisma APP

The prisma APP supports and motivates patients to adhere to treatment.

prismaTS

Therapy software for experts. For detailed data analysis and remote settings.

The most important features at a glance

- Two dynamic options for pressure adjustment in APAP mode
- Deep sleep indicator prisma RECOVER for evaluation of sleep quality
- Several convenience functions such as autoSTART, softSTART, pressure relief softPAP
- Detection of periodic breathing, RERA, snoring, hypopnea, apnea and flow limitations
- Reliable differentiation of obstructive apnea from central apnea, thanks to FOT (Forced Osciallatory Technique) ...

... and it's quiet as a whisper for the well-being and comfort of patients and partners.



prisma SOFT plus / prisma SMART plus prisma SOFT max / prisma SMART max

Accessories

prismaAQUA	WM 229490	Breathing tubes		Modules	
Humidifier, white		Breathing tube		prismaPSG	WM 29690
Set of 12 pollen filters	WM 29652	Ø 15 mm	WM 29988	PSG module	
DC Adapter 12/24 V	WM 24469	Ø 22 mm	WM 24445	prisma HUB	WM 31660
Medistrom™ Pilot-24	WM 31704	Breathing tube, autoclavak	ole	Modem for prisma CLOUD	WM 31240
lite Battery		Ø 22 mm	WM 24667	Software	
Transport Bags	Bags		VA/NA 20045	prismaTS	WM 93335
prismaBAG basic	WM 29708	HYBERNITE Superday Heated breathing tube	WM 29045	Therapy software	
prismaBAG premium	WM 29709	Ø 22 mm		complete with USB data cable	

You'll find more information about our therapy solutions, accessories and mask systems at loewensteinmedical.com

Technical data

Product class as per Directive 93/42/EEC	II a		Mean sound pressure level/Operation ISO 80601-2-70	approx. 26.5 dB(A) at 10 hPa		
Dimensions (W x H x D)	170 x 135 x 180 mm		CPAP operating pressure range	4 to 20 hPa		
Weight	ca. 1.34 kg		softSTART adjustable	5 to 45 min		
Temperature range • Operation: • Transport and Storage:	+5 °C to +40 °C -25 °C to +70 °C			4.0 hPa: 150 l/min 8.0 hPa: 174 l/min 12.0 hPa: 174 l/min 16.0 hPa: 170 l/min 20.0 hPa: 164 l/min		
Permitted humidity, operation and storage:	rel. humidity 10% to 95%, no condensation		Highest flow as per ISO 80601-2-70 Mean flow at patient connection port, 22-mm tube			
Pressure range:	700 hPa to 1060 hPa, corresponds to altitude of 3000 m above sea level automatic altitude adjustment		Stability of dynamic pressure (short-	With 22-mm breathing tube, with and without humidifier: $\Delta p \le 0.3 \text{ hPa}$ $\Delta p \le 0.4 \text{ hPa}$ $\Delta p \le 0.5 \text{ hPa}$		
Electric rating:	max. 40 VA		term accuracy)at 10 breaths/min. as per ISO 80601-2-70 at: 4.0 hPa			
System interface:	24 V DC max. 5 VA		8.0 hPa 12.0 hPa 16.0 hPa			
Power consumption in	240 V	100 V	16.0 NPa 20.0 hPa	$\Delta p \le 0.5 \text{ hPa}$ $\Delta p \le 0.6 \text{ hPa}$	Δp ≤ 0.5 hPa Δp ≤ 0.6 hPa	
 Operation (treatment): bei Bereitschaft (Standby):	0.16 A 0.035 A	0.36 A 0.061 A				
Classification as per EN 60601-1-11			Stability of static pressure (long-term accuracy) ISO 80601-2-70	Δp < 0.25 hPa +/- 3 % of measurement		
 Type of protection from electric shock: Level of protection from electric shock: Protection from damaging ingress of water and solids: 	protection class II type BF IP21		Pollen filter to 1µm to 0,3µm	Filter class E10 ≥ 99.5 % ≥ 85 %		
3			ιο σ, σμπι	≥ 00 70		

Device description Ord. Nr.

prisma SOFT plus, CPAP therapy device with Bluetooth LMT 31930 prisma SOFT max, CPAP therapy device with Bluetooth and modem LMT 31970 prisma SMART plus, APAP therapy device with Bluetooth LMT 31900 prisma SMART max, APAP therapy device with Bluetooth and modem LMT 31950

C€ 0197

Please refer to the current price list.





Löwenstein Medical Technology Kronsaalsweg 40 22525 Hamburg, Germany T. +49 40 54702-0 F. +49 40 54702-461 info@loewensteinmedical.com

Headquarters

Löwenstein Medical Arzbacher Straße 80 56130 Bad Ems, Germany T. +49 2603 9600-0 F. +49 2603 9600-50 info@loewensteinmedical.com loewensteinmedical.com







