

Full Face Mask



Performance meets comfort

Comfort is unique. Comfort is personal. Comfort is one key component to successful therapy.^{1,2} That's why we believe comfort is important across a full night's sleep, not just during the initial consultation.

F&P Vitera[™] technologies work together to provide stability, breathability, and comfort throughout the night.



90%

of patients rated the Vitera as equally or more comfortable on the face compared to their usual mask.³

Secure and stable system.


93% of trial patients rated the mask as stable, and found that the mask stayed in place during use. 80% of these patients also experienced the same or less leak compared to their usual mask.^{4,5}

Freedom to move

RollFit[™] XT, with 93% more roll than our previous generation RollFit, in combination with the stability bar, helps to keep the seal in place.⁶ The soft silicone seal rolls back and forth. Vitera's small seal is now smaller, and the large seal is now larger, than our previous generation full face mask.








Vitera headgear features VentiCool, a proprietary breathable fabric inspired by sportswear that allows more air and moisture to escape than the material used in our leading range of full face masks.*

95%

of patients rated the headgear as equally or more comfortable than their usual mask headgear.8

^{*} Internal testing shows the headgear with VentiCool technology allows 21 times more air-flow and at least 35% more moisture transfer than the headgear used in our leading full-face masks.9

Key features and benefits



Fitting guide



Step 1.

Hold the front of the mask with one hand and the headgear with the other. Ensure that one headgear clip is unhooked from the frame.



Step 2.

Place the seal onto your face and guide the headgear over your head.



Step 3.

Hook the unattached headgear clip onto the frame.



Step 4.

Gently tighten the headgear straps of the headgear starting with the bottom straps, followed by loosely fitting the blue forehead straps.

Oceania, Canada, Western Europe, Middle East (U Model)*

CODE	DESCRIPTION
VIT1SU	Vitera Mask Full Face Small
VIT1MU	Vitera Mask Full Face Medium
VIT1LU	Vitera Mask Full Face Large
VIT1SSU	Vitera Mask Small/Small (two seals)
VIT1MMU	Vitera Mask Medium/Medium (two seals)
VIT1LLU	Vitera Mask Large/Large (two seals)
VIT1MLU	Vitera Mask Medium/Large (two seals)

Asia, Eastern Europe (Z Model)*

CODE	DESCRIPTION	
VIT1SZ	Vitera Mask Full Face Small	
VIT1MZ	Vitera Mask Full Face Medium	
VIT1LZ	Vitera Mask Full Face Large	

Spare Parts

CODE	DESCRIPTION
400VIT131	Vitera Mask Full Face Mask No H/G Small*
400VIT132	Vitera Mask Full Face Mask No H/G Medium*
400VIT133	Vitera Mask Full Face Mask No H/G Large*
400VIT111	Vitera Seal Small
400VIT112	Vitera Seal Medium
400VIT113	Vitera Seal Large

CODE	DESCRIPTION
400VIT121	Vitera Headgear Small
400VIT122	Vitera Headgear Medium/Large
400FPH241	F&P Swivel
400VIT151*	Vitera Clips and Forehead Clip
400VIT171	Vitera Elbow

1. Silva, R. S. et al. An orientation session improves objective sleep quality and mask acceptance during positive airway pressure titration. Sleep Breath. (2008). doi:10.1007/s11325-007-0138-6. 2. Kryger, M., Roth, T. & Dement, W. Principles and Practice of Sleep Medicine. (Elsevier, 2010). 3. 37 out of 41 patients rated the mask as equally comfortable or more comfortable than their current mask. Internal validation study conducted on 41 participants in United States of America 2017. 4. 38 out of 41 patients rated the mask to be equally stable or more stable than their current mask. Internal validation study conducted on 41 participants in United States of America 2017. 5. 33 out of 41 patients rated the leak to be the same as or less than their current mask. Internal validation study conducted on 41 participants in United States of America 2017. 6. Internal testing shows RollFit™ XT has 93% more roll than the previous generation RollFit™. 7. 36 out of 41 patients reported they could sleep in their preferred position always or most of the time. Internal validation study conducted on 41 participants in United States of America 2017. 8. 39 out of 41 patients rated the comfort of the headgear the same or better than their current mask. Internal validation study conducted on 41 participants in United States of America 2017. 9. Internal study of F&P Vitera testing breathability of fabric with VentiCool technology against F&P Simplus™ and F&P FlexiFit™ headgear materials.

Vitera, RollFit, VentiCool, Simplus, FlexiFit and VisiBlue are trademarks of Fisher & Paykel Healthcare Limited. For patent information, see www.fphcare.com/ip



^{*}Not available in all countries.